Relaxxed

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Relaxed warmtezak

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Relaxed heat pad

EN

The heat pad can be hard and stiff upon receipt. Boil the heat pad before first use as described in this manual.

There are only three simple steps to (re)use the Relaxed heat pads.

Step 1: Click

Click the coin between your thumb and index finger. The content will react immediately and become warm.

Step 2: Use


Ready to use immediately! Place on the body part. Squeeze the heat bag for extra comfort.

Relaxxed heat pad

EN

Relaxxed heat pads give off a heat of approximately 54 °C once activated. They can be reused several hundred times.

They can be used for (among other things) warm hands / warm feet in the winter, or to relieve pain in the back or neck.

Step 3: Boiling

Boil the heat pad in a dishcloth for about 15 minutes. Wait for the water to cool.

Ready for round 2!

The heat pad is now liquid again and ready to be used again.

Caution!

1. In case of leakage: do not swallow the gel from the heat pad and avoid contact with the eyes.
2. In case of skin contact with the gel, rinse well with soap and warm water.
3. If the gel is swallowed, immediately rinse the mouth with plenty of water.
If necessary, consult a doctor.
4. In case of eye contact with the gel, rinse the eye immediately with plenty of water. (min.15 minutes).
If necessary, consult a doctor.
5. Do not leave the heat pad unattended during cooking and do not let the pan boil dry.
6. The product can be very hot after boiling!

Relaxed heat pad

7. Do not use sharp objects near the heat pocket.
8. We do not recommend using the product on areas with very sensitive skin.
9. Keep out of reach of children under 3 years old.
10. Never try to repair a damaged heat pad yourself.
11. The product and its contents are non-toxic.
12. The gel contains Sodium acetate.

Relaxxed

www.relaxxed.nl